



eve[®]

User Management & Reporting

eve[®] – The Bioprocess Platform Software

Package User Management & Reporting

eve[®] offers a clever user management system allowing you to track and document all of your bioprocess activities. You can also generate detailed reports and PowerPoint presentations from your results at the touch of a button.

User Management

- 5 predefined user roles, including a dedicated role for IT system administrators
- Computer-based access control (IP)
- Automatic logoff and password renewal
- Read-only access for guests and audits
- Highly secure passwords
- Meets technical requirements for FDA 21 CFR Part 11

Reporting

- Reports for planned batches with data ranging from device settings to batch strategy, incl. culture media and organisms
- Complete reports for ongoing and completed experiments, including audit trails, diagrams, soft sensors and more
- Reports can be exported to PDF, Word, Excel and PowerPoint at the touch of a button

Request demo and
more information at

www.infors-ht.com/eve
