












ADHÄSIVE VERSCHLUSSFOLIEN

ADHESIVE SEALING FILMS

 PERFORMANCE LINE+



REF	MAT			  	
04-095-0000	PE	80 µm	-40 °C– +104 °C	100 100 500	für / for ELISA & EIA
04-095-0050	PE	70 µm	-40 °C– +104 °C	100 100 500	ultraklar, Mikro-Kapsel-Kleber / ultra clear, non-tacky 
04-095-0060	PE	50 µm	-40 °C– +120 °C	100 100 500	ultraklar / ultra clear 
04-095-0070	PE	50 µm	-40 °C– +120 °C	100 100 500	ultraklar, für Platten mit erhöhten Rändern / ultra clear, fits raised rim plates 
04-095-0100	PE / PP	70 µm	-40 °C– +90 °C	50 50	96 klebefreie Zonen, leicht durchstechbar / 96 adhesive free zones, easy pierceable
04-096-6000	AL	60 µm	-80 °C– +120 °C	100 100 500	für Kältelagerung / for cold storage 
04-097-0000	PE	60 µm	-40 °C– +120 °C	100 100 500	für / for PCR, ELISA & EIA
04-097-5050	PE	50 µm	+4 °C– +40 °C	50 50 500	für Zellkultur (luftdurchlässig) / for tissue culture (permeable) 


04-095-0060




04-095-0070



04-095-0100



04-096-6000



04-097-0000



04-097-5050

