

# Adaptable footprint


Versatile single benches:

- ✓ Mobile with castors, load capacity up to 100 kg
- ✓ Flexibly adjustable height without castors,  
load capacity up to 300 kg
- ✓ Also available as a heavy-duty bench for point  
loads up to 200 kg or surface loads up to 400 kg



## Efficient weighing – Special balance bench

- ✓ For vibration-damped weighing
- ✓ Free-standing or for integration in a bench arrangement
- ✓ Available for sitting and standing work height





|                              | Width          |      |      |      | Height |     | Bench depth |     |
|------------------------------|----------------|------|------|------|--------|-----|-------------|-----|
|                              |                |      |      |      |        |     |             |     |
| Single bench without castors | 900            | 1200 | 1500 | 1800 | 750    | 900 | 600         | 750 |
| Single bench with castors    | 900            | 1200 | 1500 | 1800 | 750    | 900 | 600         | 750 |
| Heavy-duty bench             |                | 1200 | 1500 | 1800 | 750    | 900 | 600         | 750 |
| Balance bench, free-standing | 900            | -    | -    | -    | 750    | 900 | 600         | -   |
| Balance bench, built-in      | approx.<br>900 | -    | -    | -    | 750    | 900 | -           | -   |


Dimensions [mm]

## Flexible mounting: Bench frames and underbench cabinet frames

- ✓ A-frame variants ensure high load capacity
- ✓ C-frame variants optimise leg room
- ✓ Variants for octagonal worktops and underbench cabinets from 90° to 360°
- ✓ High load capacity with both variants:
  - C-frame: 250 kg/m<sup>2</sup> bench area
  - A-frame 350 kg/m<sup>2</sup> bench area



|         | Width |     |      |      |      | Height |     | Bench depth |     |     |
|---------|-------|-----|------|------|------|--------|-----|-------------|-----|-----|
| A-frame | 600   | 900 | 1200 | 1500 | 1800 | 720    | 870 | 600         | 750 | 900 |
| C-frame | 600   | 900 | 1200 | 1500 | 1800 | 720    | 870 | 600         | 750 | 900 |



Dimensions [mm]