

Lab furniture




INSTRUMENT TABLES QE SERIES

FLEXIBLE - MOBILE - ROBUST

The LGT-QE series

Mobile heavy-duty instrument tables with optional noise protection cabinet for vacuum pumps




Details




① Pump drawer for service



② available with single or double noise protected pump cabinet



③ Each wheel can be adjusted in height to compensate uneven floors



Details




① Monitor holder available



② Lockable heavy-duty castors



③ Optional pump status display available*



④ Pump vibration is isolated from the table frame





IMPROVE YOUR SCIENCE



SONATION GMBH
ALTE SCHULSTRASSE 39
DE-88400 BIBERACH

TEL.: +49 (0) 7351 80093-0
E-MAIL: MAIL@SONATION.COM
WWW.SONATION.COM